

# Im.PAc.T SHOULDER RETRACTOR

Model No. : R00SR27-00

Panacea's Im.PAc.T Shoulder Retractor is beneficial for shoulder clearance while treating head & neck region. Shoulder retractor plays an important role in holding the hands stiff and immobilize during the treatment. It is constructed of durable material with a large footboard and a retractor with wrist locking handles.



## Features...



Handle adjustments with hasp locking system



Indexable on a variety of couch tops



Maximum weight bearing capacity



Quick setup by fixing on indexing bar



Maximum stiffness with the retractor